

Technical Data Sheet

## Takyon™ One-Step Low Rox Probe 5X MasterMix dTTP

UF-LP5X-RT0101 • UF-LP5X-RT0501 • UF-LP5X-RT0505 • UF-LP5X-RT0510  
[0.6 mL] [5 mL] [5 x 5 mL] [10 x 5 mL]

Emerging from the combination of an optimized reaction buffer and the new efficient «Takyon™» enzyme, Takyon™ kits for Probe Assays ensure sensitivity and fast delivery of accurate and reproducible results!

### Kit contents (Table 1)

The kit UF-LP5X-RT0501 (UF-LP5X-RT0101) contains enough reagents for up to 1250 (150)- 20 µL reactions using the performant hotstart Takyon™ enzyme.

Table 1

Reagent	Volume	Description
5x MasterMix (red cap)	5 x 1 mL  0.6 mL for UF-LP5X-RT0101	Five tubes/bottles of 5x reaction mix containing e.g.: – Takyon™ DNA polymerase, – MgCl <sub>2</sub> (5.5mM final concentration), – dNTPs, – Rox Passive reference (low concentration), – Stabilizers.
50 mM MgCl <sub>2</sub> (clear cap)	1.5 ml	50 mM MgCl <sub>2</sub> solution (optional use)
Euroscript II RT/ RNase inhibitor (white cap)	250 µL	50 µL/RT 20 µL/RNase inh.
Additive (blue cap)	250 µL	Improve results on some viral and FFPE samples

### Storage conditions

5X MasterMix component of the Takyon™ One-Step Low Rox Probe 5X MasterMix dTTP should be stored between -15 °C and -25 °C and the Euroscript II RT at -65 °C to -75 °C in a constant temperature freezer. When stored under these conditions, this component is stable for 24 months (1 year for the RT). For short term storage the Takyon™ One-Step Low Rox Probe 5X MasterMix dTTP can be stored at 4 °C for 6 months.

Euroscript II RT & additive should be stored at a temperature between -65 °C and -75 °C in a constant temperature freezer. When stored under these conditions, the reagents are stable for 1 year. For short term storage the Euroscript II RT and additive can be stored between -15 °C and -25 °C for 1 month

### Procedure

- 1- Thaw all required reagents completely and put them on ice except for the Euroscript II RT and RNase inhibitor, which should be kept in the freezer until required for use. Mix all reagents well by inversion and spin them down prior to pipetting.
- 2- Prepare the reaction mix (see Table 2). Reaction set up should be done on ice. To correct for dispensing losses, prepare an excess of reaction mix (e.g. a 100-reaction mix for 96 reactions).
- 3- Add all components together, except for the template. Mix thoroughly by pipetting or inversion. Spin down.

Table 2

Component	Volume (µL)	Final Concentration
Takyon™ MasterMix	4	1x
Forward primer	2	50-900 nM <sup>1</sup>
Reverse primer	2	50-900 nM <sup>1</sup>
Probe	2	100-250 nM <sup>1</sup>
Euroscript II RT	0.2 µL	10u
Additive (optional) <sup>3</sup>	0.2 µL	
RNase-free Water	up to 17.5 µL	(volume is 20 µL minus all other components) <sup>2</sup>
Total Mix / reaction	17.5 µL <sup>2</sup>	

4- Add the reaction mix to individual reaction vials.

5- Add the template to individual reaction vials, gently mix on a magnetic stirrer and centrifuge to avoid bubbles. Negative control containing no RNA template should always be included. Optionally, a no RT-control should be set up in tubes / wells, which does not contain the EuroScript II RT/RNase Inhibitor.

6- The Takyon™ One-Step Low Rox Probe 5X MasterMix dTTP will produce consistent and sensitive results under FAST and REGULAR cycling conditions. Program the Real-Time thermocycler using the following recommended parameters (Table 3):

Table 3

	T°C	FAST cycling* FAST ramping rates! - Only on FAST cyclers	Regular Cycling Regular ramping rates!
	Time		
<b>a) Reverse transcription</b>	48°C	10 min.	10 min.
For difficult templates, increase RT step by increment of 10', up to a total of 30', to improve reaction yield			
<b>b) c-DNA amplification step:</b>			
Takyon™ activation	95 °C	3 min.	3 min.
40 Cycles			
Denaturation	95 °C	3 sec.	10 sec.
Annealing / extension	60 °C **	20 - 30 sec.	45 - 60 sec.

\* Only perform fast cycling on FAST cyclers equipped with a FAST block. Short amplicons (<120 bp) are recommended to support FAST cycling conditions. For longer amplicons or difficult templates, increase the annealing-extension time up to 40 sec.

Example of FAST cyclers: LC480, RotorGenes, ABI 7500 & 7900 with FAST block (optional), ViiA7, ABI StepOne Plus, MasterCycler ep realplex with FAST block (optional),...

\*\* The annealing temperature will vary depending on the melting temperature (T<sub>m</sub>) of the primers. Note that some FAST thermocyclers can accommodate shorter annealing steps for faster qPCR results. However some assays may require longer extension times for efficient amplification. Increase extension time by increments of 5-second, if required.

## Technical Data Sheet

### Technical information

#### Primer and probe design guidelines

##### Probes:

- Avoid runs of identical nucleotides, especially of 4 or more Gs.
- The probe T<sub>m</sub> should be 7 to 10 °C above primers T<sub>m</sub>.
- Avoid 5'-end G as it quenches the fluorophore.
- For genotyping, the position of the polymorphism should be in the centre of the probes, and the probe length should be adjusted such that each probe has the same T<sub>m</sub>.

##### Primers:

- GC content should be between 30 % and 80 % (ideally 40-60 %).
- Avoid runs of identical nucleotides, especially of 3 or more Gs or Cs at the 3' end.
- The T<sub>m</sub> should be between 58 °C and 60 °C.
- The primer should be placed as close as possible to the probe.

#### Custom assay design

The commonly used concentrations for primers and for probes are 300nM and 100nM respectively. Optimal results may require titration of primers and probes or adjustment of the primer / probe ratio. The purpose of such a process is to determine the minimum amount of primers and probe required to obtain the most sensitive results with your assay.

##### Primer titration matrix

Titrate according to the Table 4, perform qPCR and select the concentration which gives the lowest C<sub>q</sub> value. By doing this type of titration it is also possible to compensate for differences up to 2 °C in melt temperature of the primers.

Table 4: Primer titration matrix

Reverse	Forward		
	50 nM	300 nM	900 nM
50 nM	50 / 50	300 / 50	900 / 50
300 nM	50 / 300	300 / 300	900 / 300
900 nM	50 / 900	300 / 900	900 / 900

##### Primer-probe ratio matrix

Select optimal primer concentration as described in Table 4 and test with all probe concentrations described in Table 5. Select the concentration which gives the lowest C<sub>q</sub> value.

Table 5: Primer-probe ratio matrix

Opt. primers conc.	Probe		
	50 nM	100 nM	250 nM

##### MgCl<sub>2</sub> adjustment matrix

Standard MgCl<sub>2</sub> concentration is 5.5 mM but optimal MgCl<sub>2</sub> concentration can vary between assays. If necessary adjust the MgCl<sub>2</sub> concentration with the provided 50 mM MgCl<sub>2</sub> tube. Always prefer optimizing the primer and probe concentrations before the MgCl<sub>2</sub> concentration.

Adjust the amount of water if MgCl<sub>2</sub> is added to the reaction.

#### For further information please contact our Customer Help Desk:

##### For Europe:

E-mail: info@eurogentec.com

Tel: +32 4 372 76 65 • Toll-free: + 800 666 00 123

##### For USA:

E-mail: service@anaspec.com

Tel.: +1-510-791-9560 • Toll-free: +1-800-452-5530

#### FOR RESEARCH USE ONLY NOTICE TO PURCHASER: LIMITED LICENSE

Licensed dsDNA-Binding Dye Research Kits - Purchase of this product includes an immunity from suit under patents specified in the product insert to use only the amount purchased for the purchaser's own internal research. No other patent rights are conveyed expressly, by implication, or by estoppel. Further information on purchasing licenses may be obtained by contacting the Director of Licensing, Applied Biosystems, 850 Lincoln Centre Drive, Foster City, California 94404, USA.

UDG - Use of UDG employs U.S. Patents 5,035,996, 5,945,313, 5,683,896 and their foreign counterparts licensed to Eurogentec S.A. from Invitrogen Corporation.

Rox Passive Reference - Purchase of this product includes an immunity from suit to use only the amount purchased for the purchaser's own internal research. No other patent rights are conveyed expressly, by implication, or by estoppel. For information about these rights or on obtaining additional rights, please contact outlicensing@lifetech.com or Out Licensing, Life Technologies, 5791 Van Allen Way, Carlsbad, California 92008.

Takyon™ is a trademark of Eurogentec.

Note 1: Primer and probe concentrations of 300 nM & 250 nM, respectively, are recommended as starting concentrations. These concentrations will be correct for many assays, but additional optimization of the primer concentrations and primer-probe ratio may be required to obtain the best results with your primer-probe set (see table 4).

Note 2: 17.5 µL of reaction mix is added to 2.5 µL of template prior to cycling, giving a final reaction volume of 20 µL. See steps 4 and 5. These volumes, including primers & probes, can be adjusted depending on the template and reaction volume.