

EUROPE

LIEGE SCIENCE PARK • 4102 Seraing • BELGIUM • Tel.: +32 4 372 74 00 • Fax: +32 4 372 75 00
Toll-free: +800 666 00 123 • info@eurogentec.com • www.eurogentec.com

NORTH AMERICA

ANASPEC - 34801 Campus Drive • Fremont, CA 94555 • USA • Tel.: +1-510-791-9560
Toll-free: +1 800-452-5530 • Fax: +1 510-791-9572 • service@anaspec.com
www.anaspec.com

Eurogentec products are sold for research or laboratory use only and are not to be administered to humans or used for medical diagnostics.

Technical Data Sheet

Takyon™ No Rox Probe Core Kit dTTP
UF-NPCT-C0201 • UF-NPCT-C0205 • UF-NPCT-C0210
[1250 RXN - 20µL] [5 x 1250 RXN - 20µL] [10 x 1250 RXN - 20µL]

Emerging from the combination of an optimized reaction buffer and the new efficient «Takyon™» enzyme, Takyon™ kits for Probe Assays ensure sensitivity and fast delivery of accurate and reproducible results!

Storage conditions

For long term storage the Takyon™ No Rox Probe Core Kit dTTP should be stored at a temperature between -15 °C and -25 °C in a constant temperature freezer. When stored under these conditions, the components are stable for 24 months. For short term storage the Takyon™ No Rox Probe Core Kit dTTP can be stored at 4 °C for 6 months.

Kit contents (Table 1)

The kit UF-NPCT-C0201 contains enough reagents to prepare up to 1250 - 20 µL reactions using the performant hotstart Takyon™ enzyme.

Table 1

Reagent	Volume	Description
10x Buffer tube (black cap)	2 x 1.5 mL	One tube of 10x reaction buffer contains: - KCl and Tris-HCl - Stabilizers,"
50 mM MgCl ₂ (clear cap)	2 x 1.5 mL	50 mM MgCl ₂ solution (optional use)
5 mM dNTPs (green cap)	1 x 1.25 mL	A blend of dATP, dCTP, dGTP and dTTP
Takyon™ enzyme (yellow cap)	1 x 125 µL	Takyon™ enzyme (5 U/µL)
dUTP/UNG mix	1 x 330 µL	dUTP and Uracyl N-glycosylase blend for carryover prevention

Procedure

- 1- Thaw all required reagents completely and put them on ice. Mix all reagents well by inversion and spin them down prior to pipetting.
- 2- Prepare the reaction mix (see Table 2). To correct for dispensing losses, prepare an excess of reaction mix (e.g. a 100-reaction mix for 96 reactions).
- 3- Add all components together, except for the template. Mix thoroughly by pipetting or inversion. Spin down.

Table 2

Component	Volume (µL)	Final Concentration
10x reaction buffer	2	1x
Forward primer	2	50-900 nM ¹
Reverse primer	2	50-900 nM ¹
Probe	2	100-250 nM ¹
50 mM MgCl ₂	2.2	5.5 mM
5mM dNTP mix	0.8	200 µM of each dNTP
Takyon™ enzyme 5U/µL	0.1	0.02 U/µL
dUTP/UNG additive ²	0.25	Optional
Water	6.15	Water volume is 17,5 µL minus volume of all other components
Total Mix / reaction	17.5 ³	
Template or Control	2.5 ³	

- 4- Pipette either 2.5 µL of the template cDNA/DNA for your samples or 2.5 µL of the control DNA for your positive control or 2.5 µL of water/buffer for your negative control into your qPCR tubes / plate.
- 5- Add 17.5 µL of the reaction mix per well / vial, close the plate / vial and mix gently on a stirrer and spin down. Ensure that no bubbles are present in the reaction wells / vials. Reaction set up can be done at room temperature.
- 6- The Takyon™ No Rox Probe Core Kit dTTP will produce consistent and sensitive results under FAST and REGULAR cycling conditions. Program the Real-Time thermocycler using the following recommended parameters (Table 3):

Table 3

	T °C	FAST cycling* - FAST ramping rates! - Only on FAST cyclers	Regular Cycling - Regular ramping rates!
Carry over prevention optional**	50 °C**	2 min.	2 min.
Takyon™ activation	95 °C	3 min.	3 min.
40 Cycles			
Denaturation	95 °C	3 sec.	10 sec.
Annealing / extension	60 °C ***	20 - 30 sec.	45 - 60 sec.

* Only perform fast cycling on FAST cyclers equipped with a FAST block. Short amplicons (<120 bp) are recommended to support FAST cycling conditions. For longer amplicons or difficult templates, increase the annealing-extension time up to 40 sec.

Example of FAST cyclers: LC480, RotorGenes, ABI 7500 & 7900 with FAST block (optional), ViiA7, ABI StepOne Plus, MasterCycler ep realplex with FAST block (optional),...

** dUTP/UNG blend must be added to the reaction mix (see table 2).

*** The annealing temperature will vary depending on the melting temperature (T_m) of the primers. Note that some FAST thermocyclers can accommodate shorter annealing steps for faster qPCR results. However some assays may require longer extension times for efficient amplification. Increase extension time by increments of 5-second, if required.

Note 1: Primer and probe concentrations of 300 nM & 250 nM, respectively, are recommended as starting concentrations. These concentrations will be correct for many assays, but additional optimization of the primer concentrations and primer-probe ratio may be required to obtain the best results with your primer-probe set (see table 4).

Note 2: If carry over contaminations is a concern optionally add 0.25µL of the dUTP/UNG additive per 20µL of reaction.

Note 3: 17.5 µL of reaction mix is added to 2.5 µL of template/control DNA prior to cycling, giving a final reaction volume of 20 µL. See steps 4 and 5. These volumes, including primers and probes, can be adjusted depending on the template and reaction volumes.

EUROPE

LIEGE SCIENCE PARK • 4102 Seraing • BELGIUM • Tel.: +32 4 372 74 00 • Fax: +32 4 372 75 00
Toll-free: + 800 666 00 123 • info@eurogentec.com • www.eurogentec.com

NORTH AMERICA

ANASPEC - 34801 Campus Drive • Fremont, CA 94555 • USA • Tel.: +1-510-791-9560
Toll-free: +1 800-452-5530 • Fax: +1 510-791-9572 • service@anaspec.com
www.anaspec.com



Technical Data Sheet

Technical information

Primer and probe design guidelines

Probes:

- Avoid runs of identical nucleotides, especially of 4 or more Gs.
- The probe T_m should be 7 to 10 °C above primers T_m.
- Avoid 5'-end G as it quenches the fluorophore.
- For genotyping, the position of the polymorphism should be in the centre of the probes, and the probe length should be adjusted such that each probe has the same T_m.

Primers:

- GC content should be between 30 % and 80 % (ideally 40-60 %).
- Avoid runs of identical nucleotides, especially of 3 or more Gs or Cs at the 3' end.
- The T_m should be between 58 °C and 60 °C.
- The primer should be placed as close as possible to the probe.

Custom assay design

The commonly used concentrations for primers and for probes are 300 nM and 100 nM respectively. Optimal results may require titration of primers and probes or adjustment of the primer / probe ratio. The purpose of such a process is to determine the minimum amount of primers and probe required to obtain the most sensitive results with your assay.

Primer titration matrix

Titrate according to the Table 4, perform qPCR and select the concentration which gives the lowest C_q value. By doing this type of titration it is also possible to compensate for differences up to 2 °C in melt temperature of the primers.

Reverse	Forward		
	50 nM	300 nM	900 nM
50 nM	50 / 50	300 / 50	900 / 50
300 nM	50 / 300	300 / 300	900 / 300
900 nM	50 / 900	300 / 900	900 / 900

Table 4: Primer titration matrix

Primer-probe ratio matrix

Select optimal primer concentration as described in Table 4 and test with all probe concentrations described in Table 5. Select the concentration which gives the lowest C_q value.

Table 5: Primer-probe ratio matrix

Opt. primers conc.	Probe		
	50 nM	100 nM	250 nM

MgCl₂ adjustment matrix

Standard MgCl₂ concentration is 5.5 mM but optimal MgCl₂ concentration can vary between assays. If necessary adjust the MgCl₂ concentration with the provided 50 mM MgCl₂ tube. Always prefer optimizing the primer and probe concentrations before the MgCl₂ concentration.
Adjust the amount of water if MgCl₂ is added to the reaction.

For further information please contact our Customer Help Desk:

For Europe:

E-mail: info@eurogentec.com
Tel: +32 4 372 76 65 • Toll-free: + 800 666 00 123

For USA:

E-mail: service@anaspec.com
Tel.: +1-510-791-9560 • Toll-free: +1-800-452-5530

FOR RESEARCH USE ONLY NOTICE TO PURCHASER: LIMITED LICENSE

Licensed dsDNA-Binding Dye Research Kits - Purchase of this product includes an immunity from suit under patents specified in the product insert to use only the amount purchased for the purchaser's own internal research. No other patent rights are conveyed expressly, by implication, or by estoppel. Further information on purchasing licenses may be obtained by contacting the Director of Licensing, Applied Biosystems, 850 Lincoln Centre Drive, Foster City, California 94404, USA.

UDG - Use of UDG employs U.S. Patents 5,035,996, 5,945,313, 5,683,896 and their foreign counterparts licensed to Eurogentec S.A. from Invitrogen Corporation.

Rox Passive Reference - Purchase of this product includes an immunity from suit to use only the amount purchased for the purchaser's own internal research. No other patent rights are conveyed expressly, by implication, or by estoppel. For information about these rights or on obtaining additional rights, please contact outlicensing@lifetech.com or Out Licensing, Life Technologies, 5791 Van Allen Way, Carlsbad, California 92008.

Takyon™ is a trademark of Eurogentec.

SYBR® Green I nucleic acid gel stain. This product is provided under a technology license from Life Technologies Corporation to Eurogentec SA. The sale of this product is conditioned on the buyer using the product only in research conducted by the buyer, and not to sell or otherwise commercialize this product or its components (1) in manufacturing; (2) to provide a service, information, or data in exchange for payment; (3) for therapeutic, diagnostic or prophylactic purposes; or (4) for resale, not even if resold for use in research. For information on purchasing a license to this product for purposes other than research, contact Life Technologies Corporation, 5791 Van Allen Way, Carlsbad, CA 92008 USA or outlicensing@lifetech.com.