

OLIGO NAME	CATALOG #	CHEMISTRY	BASES	DATE	
HL647-T10	UN-CT030-005	DNA	10	03-Jul-18	
SEQUENCE (5' → 3')					
5' HiLyte™ Fluor 647 - TTT-TTT-TTT-T 3'					
BACKBONE (BASES + LINKAGES)					
A	0	GC (%) ¹	0.0	Quantity (OD@260nm)	0.50
C	0	MW (g.mol ⁻¹) ¹	4372.7	(nmol)	5.0
G	0	Ext. coeff. (l.mol ⁻¹ .cm ⁻¹) ¹	99300	(µg)	21.9
T	10	Tm (°C) ¹	20.0		
Others	1				
¹ Theoretical values					
PURIFICATION(S)					
HPLC-RP					
QUALITY CONTROL(S)					
MS					Passed

DELIVERY CONDITION
Room Temperature
RECONSTITUTION
To make up to 100 µM concentration, add 50 µl of recommended buffer. To make up to 20 µM concentration, add 250 µl of recommended buffer. To make up to 5 µM concentration, add 1000 µl of recommended buffer.
LONG TERM STORAGE
Store at -20 °C

COMMENTS		
General: Tm % GC calculated with 50 mM NaCl and 0 % formamide (PCR conditions) Specific: Stable until 03-Jan-20 under dry condition at room temperature, QC available on request		
TECHNICAL SUPPORT		
If you have any questions feel free to call our Oligo Centre		
<table border="0" style="width: 100%;"> <tr> <td style="width: 50%; vertical-align: top;"> EUROPE ☎ 00 800 666 00 123 (European toll free number), ✉ oligo@eurogentec.com Kaneka Eurogentec S.A. Liège Science Park Rue Bois Saint-Jean 5 - 4102 SERAINING BELGIUM Tel: +32(0)4 372 74 00 - Fax: +32(0)4 264 07 88 E-mail: info@eurogentec.com Web: www.eurogentec.com RPM Liège T.V.A.-(BE)-0427.348.346 - ING Belgique Bank - IBAN: BE86 3400 2118 6050 BIC: BBRUBEBB </td> <td style="width: 50%; vertical-align: top;"> NORTH AMERICA ☎ +1 (510) 896-1350 ✉ oligo.na@egt-biotech.com EGT North America, Inc. 34801 Campus Drive - Fremont, CA 94555 Tel.: +1 (510) 896 -1350 * Fax: +1 (510) 791-9572 E-mail: ivd.na@egt-biotech.com * Web: www.egt-biotech.com </td> </tr> </table>	EUROPE ☎ 00 800 666 00 123 (European toll free number), ✉ oligo@eurogentec.com Kaneka Eurogentec S.A. Liège Science Park Rue Bois Saint-Jean 5 - 4102 SERAINING BELGIUM Tel: +32(0)4 372 74 00 - Fax: +32(0)4 264 07 88 E-mail: info@eurogentec.com Web: www.eurogentec.com RPM Liège T.V.A.-(BE)-0427.348.346 - ING Belgique Bank - IBAN: BE86 3400 2118 6050 BIC: BBRUBEBB	NORTH AMERICA ☎ +1 (510) 896-1350 ✉ oligo.na@egt-biotech.com EGT North America, Inc. 34801 Campus Drive - Fremont, CA 94555 Tel.: +1 (510) 896 -1350 * Fax: +1 (510) 791-9572 E-mail: ivd.na@egt-biotech.com * Web: www.egt-biotech.com
EUROPE ☎ 00 800 666 00 123 (European toll free number), ✉ oligo@eurogentec.com Kaneka Eurogentec S.A. Liège Science Park Rue Bois Saint-Jean 5 - 4102 SERAINING BELGIUM Tel: +32(0)4 372 74 00 - Fax: +32(0)4 264 07 88 E-mail: info@eurogentec.com Web: www.eurogentec.com RPM Liège T.V.A.-(BE)-0427.348.346 - ING Belgique Bank - IBAN: BE86 3400 2118 6050 BIC: BBRUBEBB	NORTH AMERICA ☎ +1 (510) 896-1350 ✉ oligo.na@egt-biotech.com EGT North America, Inc. 34801 Campus Drive - Fremont, CA 94555 Tel.: +1 (510) 896 -1350 * Fax: +1 (510) 791-9572 E-mail: ivd.na@egt-biotech.com * Web: www.egt-biotech.com	