

EUROPE

LIEGE SCIENCE PARK • 4102 Seraing • BELGIUM • Tel.: +32 4 372 74 00
• Fax: +32 4 372 75 00 • Toll-free: + 800 666 00 123 • info@eurogentec.com
• www.eurogentec.com



NORTH AMERICA

ANASPEC - 34801 Campus Drive • Fremont, CA 94555 • USA • Tel.: +1-510-791-9560
• Toll-free: +1 800-452-5530 • Fax: +1 510-791-9572 • service@anaspec.com
• www.anaspec.com

Technical Data Sheet

Electroporation Cuvettes CE-0001-50 • CE-0002-50 • CE-0002-50S • CE-0004-50

Eurogentec products are sold for research or laboratory use only and not for medical diagnostics.

Description

Eurogentec electroporation cuvettes are moulded in high quality polycarbonate, with extremely precise and uniform electrode gaps.

The electroporation cuvettes are compatible with most existing electroporation systems and are available in 3 sizes:

- 1 mm (black cap),
- 2 mm (yellow cap)
- 4 mm (blue cap)

Product Specification

Material Polished Aluminium

Cat #	Cap color	Gap (mm)	Min. Volume (µL)	Max. Volume (µL)
CE-0001-50	Black	1	10	90
CE-0002-50	Yellow	2	20	400
CE-0002-50S	Yellow	2	20	400
CE-0004-50	Blue	4	80	800

Sterility Gamma irradiated sterile
Gap tolerance +/- 2%
Max electricity field 25kV/cm
Packaging 50, individually sealed

Storage

The electroporation cuvettes can be stored at room temperature.

Shipping

The electroporation cuvettes are shipped at room temperature.

For further information please contact our Customer Help Desk:

For Europe:
E-mail: info@eurogentec.com
Tel: +32 4 372 76 65

For USA:
E-mail: service@anaspec.com
Tel: +1 510 791 9572