

EUROPE

LIEGE SCIENCE PARK • 4102 Seraing • BELGIUM • Tel.: +32 4 372 74 00
• Fax: +32 4 372 75 00 • Toll-free: + 800 666 00 123 • info@eurogentec.com
• www.eurogentec.com



NORTH AMERICA

ANASPEC - 34801 Campus Drive • Fremont, CA 94555 • USA • Tel.: +1-510-791-9560
• Toll-free: +1 800-452-5530 • Fax: +1 510-791-9572 • service@anaspec.com
• www.anaspec.com

Technical Data Sheet

iD Sample Buffer 5X ID-WCSAB1-005

Eurogentec products are sold for research or laboratory use only and are not to be administered to humans or used for medical diagnostics.

Product

Cat #	Description	Size
ID-WCSAB1-005	iD Sample Buffer 5X	5 mL

Description

iD Sample Buffer 5X is used as a tracking dye for SDS PAGE gel loading. This buffer contains SDS and is suitable for denaturing gel electrophoresis.

Storage

Store at -20°C

Application

Dilute 1:4 with sample before loading.

Note: 5X Sample Buffer should be brought to room temperature prior to use.

Mix well before loading gel.

For further information please contact our Customer Help Desk:

For Europe:

E-mail: info@eurogentec.com
Tel: +32 4 372 76 65

For USA:

E-mail: service@anaspec.com
Tel: +1 510 791 9560